



LOGIQ™ P9

Performance Series

Make it easy. Make it your own.



Imagination at work

LOGIQ is a trademark of General Electric Company or one of its subsidiaries.

Global
January 2015
JB26907XE



Challenging ultrasound environment

Time and quality challenges

- Growing patient volume
- Need for fast, accurate diagnosis

Operator essentials

- User-friendly equipment
- Avoid work related injuries

Economic pressures

- Budget constraints
- Declining reimbursement
- ROI and life cycle costs

Physical environment

- Cramped exam rooms
- Increasing need for portable exams

Patient expectations

- Well-informed patients
- Comfort and safety





GE Ultrasound innovation

Clinical focus

LOGIQ™

General
Imaging



Voluson™

Women's
Health



Vivid™

Cardio-
vascular



Venue™

Point of Care



Vscan™

Primary Care



ViewPoint™

Ultrasound IT



**Invenia™
ABUS**

Breast
Screening



**Extraordinary
images**

Easy workflow
Image fusion
and scan
automation

**Extraordinary
images**

Exceptional 4D
Advanced
automation

**Extraordinary
images**

4D heart in a
single beat
Advanced
quantitative tools

**Extraordinary
images**

Easy to use
Easy to clean

**Extraordinary
images**

“Take a look”
Size/weight of a
soda can

**Ultrasound-specific
productivity**

Easy Workflow
Plug-ins for popular
PACS

**Designed for
screening**

Advanced
automation for
fast-paced
workflow





GE Ultrasound innovation

General Imaging products

Compact

Venue™ Series LOGIQ™ e Vscan™



Invenia™
ABUS



LOGIQ E9
with XDclear™



LOGIQ P7/LOGIQ
P9



LOGIQ S7



LOGIQ S8





LOGIQ™ P9

Make it easy. Make it your own.

Easy to use systems with fast performance and personalized workflow



Personalized

Customizable settings and workflow that help improve your ease of use



Patient-Centric

Solutions that enable improved focus on patients and help enhance your confidence and care delivery



Practical

Smart ergonomics and solutions to help improve your productivity and investment value



Personalized



LOGIQ™ P9

Make it easy. Make it your own.



My Page for customized personal preferences

- Personalize workflow by preset
- Configurable measurement, comments, body patterns, and functions
- Individual user logon

Extended customization for enhanced flexibility

- User defined keys in control panels and touch panels
- Configurable mode keys
- Configurable smart keys
- Digital TGC preset

Simplified operator panel for fast exams

- 80% reduction of physical buttons and controls*
- 10.4 inch touch panel
- Digital TGC and keyboard



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.
*Internal GE engineering study using standardized protocols for an abdominal exam compared with prior version GE LOGIQ P6 ultrasound system.



Global
January 2015
JB26907XE

Patient-Centric





LOGIQ™ P9

Helps you improve patient care

Achieve high image quality with ease

- Migration of [Agile Acoustic Architecture](#) from our leadership system into a lightweight portable design
- Wide assortment of [LOGIQ advanced probes](#) for high performance in a wide range of applications
- Features to help enhance diagnostic information such as Speckle Reduction Imaging, [CrossXBeam™](#) and [B-steer+²](#)

Advanced tools² to apply on difficult cases

- Vascular features as [Contrast Imaging¹](#) and [B-Flow™](#),
- [Elastography](#)
- [LOGIQ View](#)
- [3D/4D Ultrasound](#)
- [Anatomical M-Mode](#)

Addressing a variety of clinical needs

- Comprehensive solutions to address your needs in a wide range of clinical applications: Abdominal, Small Parts, Breast, Musculoskeletal, Vascular, Cardiac, OB/GYN



LOGIQ, CrossXBeam and B-Flow are trademarks of General Electric Company or one of its subsidiaries.

1. CEUS is not approved or cleared by the U.S. FDA.
2. Refer to the LOGIQ P9 brochure for a list of purchasable options



Global
January 2015
JB26907XE


Radiology



Image uniformity with C1-5-RS from near to far field



Thyroid lesion vascularization using Power Doppler imaging and ML6-15-RS



Internal carotid artery plaque detailed representation using 9L-RS




Abdominal



Uniform liver imaging with C1-5-RS from near to far field




Vascularization of the spleen using B-Flow™ and 9L-RS




Renal artery hemodynamics with PW using C1-5-RS




Vascular



Clear visualization of the IMT of the CCA using 9L-RS




True hemodynamics seen in vascular imaging with B-Flow™ Color and 9L-RS




ICA hemodynamics with PW using 9L-RS




Small Parts



Lesion differentiation in breast using ML6-15-RS




Stiffness representation of breast lesion using Elastography option with ML6-15-RS




Thyroid lesion vascularization using B-Flow™ Color and 12L-RS



Musculoskeletal



Representation of the subscapularis tendon using ML6-15-RS




Biceps tendon imaging with ML6-15-RS




Superficial imaging of MTP joint with ML6-15-RS




Cardiology



Aortic valve regurgitation flow in apical long axis view using 3Sc-RS



Mitral valve M-Mode imaging with 3Sc-RS




Assessment of wall motion using TVI/TVD option with 3Sc-RS



OB/GYN



3D/4D imaging of the fetus with RAB2-6-RS



Endocavitary imaging of the uterus using E8C-RS



Image uniformity with RAB2-6-RS from near to far field



Practical





LOGIQ™ P9 helps you...

Ergonomics on demand

An innovative ergonomic solution addressing your request for a slim, fully adjustable system to help increase your scanning comfort

Solutions¹ to help you improve exam consistency and productivity

- [Auto Image Optimization \(ATO, ASO\)](#)
- Reset key
- Breast and Thyroid [productivity packages](#)
- [Compare Assistant](#)
- [Scan Assistant](#)
- [GE Raw Data](#)
- [Auto IMT](#)
- [AutoEF](#)
- [Multi Modality Q/R](#)
- [Measure Assistant](#)



Extended portability to address your mobility needs

- [Power Assistant](#) and wireless capability



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.
1. Refer to the LOGIQ P9 brochure for a list of purchasable options



Global
January 2015
JB26907XE

Customer care



LOGIQ™ P9

Make it easy. Make it your own.


Tailored solutions¹ to help you address your daily challenges

Education

- My Trainer – on-board training tools
- LOGIQ™ Club community and product educational services
- Access to additional purchasable online and in-person education

Service and Care

- Three years of service included with purchase
- Global remote support (InSite™ Exc)
- Dedicated financial solutions
- Probe care
- www.gehealthcare.com/transducers





LOGIQ™ P9

Performance Series

Make it easy. Make it your own.

Its sleek, lightweight design incorporates powerful technologies while enhancing workflow and user experience.

All at an affordable price that helps meet budget expectations.



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.

Global
January 2015
JB26907XE






My Page

Configurable on each My Page

- GE innovative personalization feature
- Preset depended (Abd/Carotid/...has My Page)
- Each use case (preset) has My Page

Abd

Carotid



OB

Breast

Measurements
Comments
Body patterns
Functions



Additional configurable capability



1 5 physical user defined keys

2 Freeze & P1 Exchangeable

5 digital user defined keys



5 configurable mode keys

Configurable smart keys

Imaging Mode	Function
Live B, B-Flow™	Focal Zone Up/Down
	Frequency Up/Down
CF, TVI, B-Flow Color, PDI	Box Steer
	Scale (PRF) Up/Down
PWD, CW, TVD	Baseline Up/Down
	Scale (PRF) Up/Down
M/D cursor	Sample volume size

3



Traditional Ultrasound

Rigid assumptions about how sounds interacts with the body

- Speed of sound is constant (1540 m/s)
- Sound attenuation is constant
- Ultrasound beam is pencil shaped
- Vessels are straight
- Blood flow is laminar



Looking at the body
as a phantom



Reality is not so simple:

- Sound attenuation and velocity vary
- Ultrasound beam is a volume
- Vessels are tortuous
- Blood flow can be turbulent

Simplifying assumptions can lead to:

- Poor image quality
- Distortion of image geometry
- Lots of adjustments to optimize images





Agile Acoustic Architecture

- Flexible clinically-based mathematical models of the body
 - Dynamically optimizes image acquisition for every body type
 - Including the reality of body types
-
- Speed of sound is variable: 1450-1560 m/s
 - Sound attenuation varies based on tissue type
 - Ultrasound beam is dynamic

Agile Acoustic Architecture

Powerful Acquisition System

Dynamic Mathematical Model



P-Agile Acoustic Architecture

Agile Architecture



Leadership core technology for advanced clinical applications

Migration of technology

P-Agile Architecture

Compact size and technology to fit P-series customers workflow

Agile System



P-Agile System






Agile Acoustic Architecture

Helping improve your diagnostic confidence

Clinical benefits

- Image uniformity
- Spatial resolution
- High frequency imaging at depth
- Reduced keystrokes to help optimize the image




LOGIQ™ P-Series Probes

Innovative technologies



Conventional Array

Active Matrix Array



ML6-15-RS

High frequency linear probe

Matrix arrays provide multiple rows of crystals, allow focusing in the near, mid and far field.

Great spatial resolution and image uniformity from near to far field

Footprint 50 mm

Ultra-high Doppler & Color frequency for excellent slow flow sensitivity




LOGIQ™ P-Series Probes


Acoustic Amplifier



Traditional Technology



Acoustic Amplifier Technology



Recaptures the unused energy that passes through the crystal

Improves sensitivity axial resolution & penetration over traditional probes



LOGIQ™ P9 Probes

Wide range of applications




Clinical application	LOGIQ P9 Probes
Abdominal	C1-5-RS, 9L-RS, 3Sc-RS, RAB2-6-RS
Vascular	ML6-15-RS, 9L-RS, 12L-RS, L8-18i-RS, P8D
Cardiology	3Sc-RS, 6S-RS
OB/GYN	RAB2-6-RS, C1-5-RS, E8C-RS
Breast	ML6-15-RS, 12L-RS
Small Parts	ML6-15-RS, 12L-RS, L8-18i-RS
Pediatrics/Neonatal	8C-RS, 3Sc-RS, 6S-RS, 9L-RS, ML6-15-RS
Urology	C1-5-RS, E8C-RS, RAB2-6-RS
Intraoperative	L8-18i-RS



CrossXBeam™

Spatial compounding imaging

- Provides 3, 5, 7 of spatial compounding
- Live side-by-side dual view display
- Compatible with:
 - Color Mode
 - PW
 - SRI-HD
 - Coded Harmonic Imaging
 - Virtual Convex



Help increase clinical confidence in all imaging modes




B-Steer+


B-Steer+ ¹ enables enhanced visualization of the needles structure during interventional procedures, helping improve user confidence and exam accuracy.

Highlights:

- Up to 12 selectable steering angles available (six each direction)
- Separate gain control for needle reflection
- Available on all linear probes
- Quick one-button operation



Virtual Convex



- Provides a convex field of view
- Compatible with CrossXBeam™
- Available on linear and phased array probes

View large anatomy in greater detail



Contrast Imaging^{1,2}

Amplitude Modulation Technique



- Great penetration and contrast sensitivity
- Excellent tissue suppression
- High image uniformity

Imaging Modes

- Dual or single display
- Hybrid contrast
- Accumulation

Features

- Dual caliper, dual timer
- TIC Q-analysis package
- Retrospective/prospective storage
- One button background storage




Available for C1-5-RS (Abdominal studies)

Contrast Imaging^{1,2}



Time Intensity Curve (TIC)

- Raw data processing for Contrast uptake
- Q-Analysis of both compressed and uncompressed data
- Up to 8 selectable ROI's
- Up to 10 parameters
- Ellipsoid or manual ROI tracing
- Anchor tracking function
- Automated motion tracking
- Automatic enabling/disabling of frames
- Trace export in ASCII format



TIC/Q-Analysis



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.

1. CEUS is not approved or cleared by the U.S. FDA.
2. Refer to the LOGIQ P9 brochure for a list of purchasable options.



Global
January 2015
JB26907XE


B-Flow™/B-Flow Color

A GE innovation¹

- Does not use Doppler processing
- Based on GE Patented coded technology
- Display real hemodynamics
- Direct visualization of blood reflectors

True hemodynamics and anatomy

- Dynamic appearance of flow
- Minimal tissue overwrite compared to Color Flow
- Excellent control of flash artifacts compared to Color Flow
- Clearly visible background image
- Independent mode control



B-Flow™/B-Flow Color


Benefits over Doppler

- No tissue overwrite
- No impact on frame rate
- Less angle dependency
- High resolution Imaging
- Background image clearly visible

Available on probes

LOGIQ P9

- 9L-RS
- 12L-RS
- ML6-15-RS
- L8-18i-RS



Strain Elastography¹

Strain imaging technology requiring a light manual compression or patient breath to perform tissue deformation. A qualitative and Semi-Quantitative² solution.


Highlights:

- High sensitivity and persistence
- Consistent pattern
- User selectable color maps
- Dual measurements
- User support by pressure quality bar and graph


Available on the following probes:

LOGIQ P9

- C1-5-RS
- ML6-15-RS
- 9L-RS
- 12L-RS



This soft lesion deforms under pressure.




This hard lesion does not deform under pressure.

$$\text{Elasticity} = \frac{\text{Stress}}{\text{Strain}}$$

Stress = axial force applied to lesion

Strain = tissue deformation due to applied stress



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.

1. Refer to the LOGIQ P9 brochure for a list of purchasable options

2. Not cleared by the US FDA and not for sale in the USA



Strain Elastography¹

Semi-Quantification²

E-Index



- Selected ROI's elasticity value
- Round or manual tracing of ROI
- E-index range from 0 – 6
- Based on GE Raw Data processing

E-Ratio

- Calculated E-Index between ROI's (up to 8)
- Representing relative stiffness

Q-Analysis over multi-frame acquisition

- Automatic skip of low quality frames
- Anchor function
- Trace export in ASCII format




LOGIQ View



LOGIQ™ View¹ allows a real-time representation of long anatomical areas. (e.g. Achilles' Tendon)

- Available on all probes
- Combine with CrossXBeam™ at linear probe
- Auto detection of scan direction
- Up to 60 cm scan length



1. Refer to the LOGIQ P9 brochure for a list of purchasable options
LOGIQ and CrossXBeam are trademarks of General Electric Company or one of its subsidiaries.



Global
January 2015
JB26907XE

3D/4D Ultrasound¹

GE Volume Imaging with 3D/4D dedicated abdominal probe

Highlights:

- Easy, quick, reproducible
- High volume data accuracy
- Comprehensive settings

Volume Modes:

- Multi-planar Imaging
- Surface rendering
- [TUI – Tomographic Ultrasound Imaging](#)
- [VCI – Volume Contrast Imaging](#)
- [Vocal- Volume Calculation](#)

Probes:

- RAB2-6-RS (Abdominal)



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.

1. Refer to the LOGIQ P9 brochure for a list of purchasable options



Global
January 2015
JB26907XE


3D/4D Ultrasound¹

TUI-Tomographic Ultrasound Imaging

Visualization mode that presents data as parallel slices (planes) through acquired Volume dataset. It works with CFM/PDI and SRI-HD.

Highlights:

- Information consistent to CT & MR format
- 3D static with Colour
- Up to 9 slices, with user selectable distance (min 0.5 mm, step by 0.1 mm) and angle
- Top left held as reference image
- Works with SRI




3D/4D Ultrasound¹

VCI – Volume Contrast Imaging

VCI is a volume acquisition technique enhancing B-Mode contrast resolution and speckle suppression.

Highlights:

Help improve assessment of lesions size, margins and internal structures for comprehensive patient management.



3D/4D Ultrasound¹

Vocal

Vocal is a 'Volume computer-aided analysis' based on a volume acquisition enabling fast and accurate volume calculations.

Highlights:

- Manual, semi-automatic or automatic borders definition
- Basic measurements as length, angle and area
- Easy corrections and contour modifications
- Can be used with any lesion or volume to measure



Flow Quantification¹

2D CFM/PDI Quantitative assessment of vascular feeding in a selected ROI

Highlights:

- Consistent, repeatable and objective measurement
- Can help in treatment planning and monitoring protocols
- Provides data to support outcome measurements

Features:

- Up to 8 selectable ROI's
- Analysis over 4/5 heart beat cycles
- Automatic or manual ROI tracing
- "Save ROI" feature for monitoring
- Manual disabling & enabling of frames
- Export traces in ASCII format



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.


1. Refer to the LOGIQ P9 brochure for a list of purchasable options



Global
January 2015
JB26907XE

AMM – Anatomical M-Mode¹

- M-Mode cursor adjustable at any plane to be vertical to myocardium
- Compatible with
 - Live image
 - Stored image – Raw data
 - Color Flow Mode



Evaluate heart function with additional data



Ergonomic Adaptable System¹

USB ports (monitor side)



USB ports (on the top of touch panel)



21.5" LCD monitor

10.4" LCD Touch Panel



Rear handle



Gel warmer



4 active probe ports

Small physical footprint



1. Refer to the LOGIQ P9 brochure for a list of purchasable options



Ergonomic Adaptable System¹

Peripheral options



Multi-purpose holder
(L & R reversible)

Physical
A/N keyboard



Articulating
Arm

Rear Tray

Paper Tray

High Cabinet

Drawer

Side Tray (L & R)



1. Refer to the LOGIQ P9 brochure for a list of purchasable options



Ergonomic Adaptable System




Rotate
+/- 30
degrees


Up/Down
810 - 910
mm



Ergonomic Adaptable System



Ergonomic Adaptable System



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.



Global
January 2015
JB26907XE

Intuitive interface

- LOGIQ™ design operator panel
- 10.4" touch screen with simplified User Interface
- User-friendly keys and backlight design



System adjustability



LCD monitor with articulation arm:




System adjustability

Up and down:



Compact



Connectivity



Output Ports

- USB
- Composite
- HDMI
- S-Video
- Ethernet

Network storage

- DICOM®
- SaveAs
- MPEGView
- Report Save As




Auto Optimization

One Button Press




Auto TGC

Continuous 2D auto TGC control




ATO (Auto Tissue Optimization)

Optimize B-Mode image to help improve contrast resolution



ASO (Auto Spectral Optimization)

Baseline and Scale (PRF)



Scan Assistant¹

- The personal assistant to user's exam
- Up to 63% time reduction²
- Up to 87% keystroke reduction²
- Help increase exam consistency
- Compatible with LOGIQ™ E9, LOGIQ S8 and LOGIQ S7 workflows

Feature Highlights


Initiates and completes user selected required measurements

Automatically steers Color Doppler

Automatically sets up Imaging controls and modes

Automatically inserts comments

Automatically reorders images to reader's preference



1. Refer to the LOGIQ P9 brochure for a list of purchasable options
 2. Internal GE engineering study using standardized protocols for an abdominal exam compared with prior version GE LOGIQ P6 ultrasound system.
 LOGIQ is a trademark of General Electric Company or one of its subsidiaries



Compare Assistant¹

Help streamline comparison to prior exams

Opportunity

Many ultrasound exams are follow ups to prior ultrasound or other modality exams.

Goal

Drive productivity for acquiring and reading the exam by designing a workflow that uses prior exam data.

Result

A quick image comparison or a replicated prior exam.



1. Refer to the LOGIQ P9 brochure for a list of purchasable options



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.

1. Refer to the LOGIQ P9 brochure for a list of purchasable options



Global
January 2015
JB26907XE

Compare Assistant¹

Help streamline comparison to prior exams

At the scanner...

Easy access to past exam data on the scanner



Side-by-side compare and store past exam image to today's image




Set image setting of live scanning to match past exam image²



Create entire new exam to match old exam

At the reading station...



Help reduce the time spent to find, open, sort, compare to prior exams



1. Refer to the LOGIQ P9 brochure for a list of purchasable options
2. Automated when LOGIQ P9 data is used



Measure Assistant¹

Designed to work in OB exams

In Scan Assistant:

- Freeze on anatomy
- Measurement auto applied
- Print to accept or easily edit if needed



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.

1. Refer to the LOGIQ P9 brochure for a list of purchasable options




Global
January 2015
JB26907XE

Measure Assistant¹

Designed to work in breast exams

Few simple steps to assist the user with breast measurements

- User bounds lesion with ROI
- System auto traces, generates height and width
- User prints to accept or easily edits as needed
- Messages on the status bar help guide the user



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.

1. Refer to the LOGIQ P9 brochure for a list of purchasable options



Global
January 2015
JB26907XE

Breast Productivity¹


Measurement package

A dedicated breast-specific measurement package that allows user to:

- Make labeling, measuring and describing lesion easy
- Leverage the BI-RADS[®] lexicon criteria/assessment
- Organizes multiple measurements into a convenient worksheet
- Send results via DICOM[®] SR



Directly from BI-RADS lexicon




Thyroid Productivity¹


Measurement package

Thyroid-specific measurement package

- Enables labeling, measuring and describing nodules, lymph nodes and parathyroids
- Multiple measurements can be organized into a convenient worksheet and sends results via DICOM® SR



Show features




Worksheet



Auto IMT¹

Auto IMT is an automated method of measuring the intima media thickness of the CCA or ICA from multiple samples across a user defined length.

- Simple and easy to operate
- Direct export of measurements to a worksheet and report page
- Including ECG trigger to help increase consistency
- Save offset distance and IMT measurement lengths to help increase reproducibility



An efficient, reproducible method of carotid artery analysis



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.

1. Refer to the LOGIQ P9 brochure for a list of purchasable options





Global
January 2015
JB26907XE

AutoEF¹



Automated Ejection Fraction (AutoEF) is a semi-automatic measurement tool used for measurement of the global EF (Ejection fraction).

- The AutoEF tool tracks and calculates the myocardial tissue deformation based on feature tracking on
- B-Mode cine loops
- AutoEF is performed on either one or both apical
- 4-chamber or 2-chamber views, in any order
- Result is presented as Ejection Fraction value for each view and average Ejection Fraction for the whole LV. All values are stored to the worksheet after the results are approved



LOGIQ is a trademark of General Electric Company or one of its subsidiaries.

1. Refer to the LOGIQ P9 brochure for a list of purchasable options



Global
January 2015
JB26907XE

Power Assistant¹

Power Assistant is an innovative solution that provides the system battery power during transport to help decrease system shut-down and reboot time – helping achieve excellent productivity for excellent portable exams.

Highlights

- Prompt In & out of battery operation mode (lasting up to 20 min.) to help improve system's portability
- Simple plug in and out operation
- System safely shuts down automatically before battery runs out
- Wireless LAN capability

Notes

- Scanning not available when in battery operation mode
- Battery charges only when system power is on



1. Refer to the LOGIQ P9 brochure for a list of purchasable options



Raw Data

The foundation for simplified workflow

Raw Data capture enables to build a thorough exam while helping reduce scan time. This **proprietary Raw Data format** from GE Healthcare captures data earlier in the image processing chain enabling users to make changes to the data during or even after the exam has ended.

Room too bright?
Adjust gain later

Difficult vascular patient?
Adjust baseline shift,
sweep speed later

Forgot annotations?
Easily add them later

Delicate NICU patient?
Acquire quickly then
virtually rescan later





Raw Data


TruScan™ Architecture

Capturing Raw Data early in the image chain

Typical System



LOGIQ



Raw Data processing with TruScan architecture enabling “Virtual Rescanning”




Raw Data processing

Original Acoustic Data are stored before Scan Converting in a GE “Raw” format to be easily accessed and re-processed any time after the exam completion.

Highlights:


- Sub-optimal studies can be optimized
- Measurements can be re-done and reports regenerated
- All Imaging control parameters can be changed as:
 - B-Mode: Gain, DR, AO, Zoom, SRI...
 - CFM: Gain, Threshold, DualView...
 - PW: Baseline, Invert, Angle, Gain...



Multi-modality Query Retrieve¹

Live Ultrasound comparison with retrieved CT or MRI Volume dataset

- Potential time saver
- Simplified workflow
- Complementarity of information
 - Volume data set retrieval
 - Select desired image plane by scrolling
- Helps improve diagnostic confidence



1. Refer to the LOGIQ P9 brochure for a list of purchasable options
LOGIQ is a trademark of General Electric Company or one of its subsidiaries.





My Trainer – on-board training modules

The image shows three parts of the My Trainer interface:

- Left Panel:** A vertical menu with a "Bookmark" button at the top. The menu items are: Getting started, Probe management, Ergonomics, Monitor arm movement, Console movement, Brakes-footrest, System transportation, Utility basics, Connectivity, Use your system, and FAQs. A large blue arrow labeled "Expandable menus" points to this menu. At the bottom, there is a "Contents" section with a list: "Get started", "Use your system", and "FAQs".
- Middle Panel:** A "Welcome to My Trainer! For LOGIQ9" screen. It features a "Bookmark" button and a smaller version of the menu from the left panel.
- Right Panel:** A training module titled "Monitor arm movement (1/2)". It includes a video of a person adjusting a monitor arm. A blue box labeled "Contents Pictures or Videos" points to the video. Below the video is a "Description" section. A blue box labeled "Contents select buttons (four): First, Previous, Next, Last" points to a set of four navigation buttons below the video.

Description
Articulating arm: Turn the release knob counter clockwise to unlock the LCD monitor. The LCD monitor can be moved freely in all directions. Grab the bottom of the LCD monitor when you adjust the position of the monitor and monitor arm.





LOGIQ™ Club users community

www.LOGIQClub.net

Access to the Club websites

- Local content and events
- Clinical in-person courses
- Application tips & tricks
- Publications, cases & papers
- Clinical cases & technical presentations
- Downloadable product educational materials and DVDs

Personalized Mailings and Newsletters

Learn about new ultrasound products and software upgrades

User Days and Lounges

- Learn best practices from specialists around the globe
- Discuss and exchange information with ultrasound users worldwide





Service. Designed for peace of mind.

With three years of coverage¹, the LOGIQ™ P9 system helps provide you value and peace of mind right from the start.

LOGIQ P9 system is also empowered by InSite™* with Express Connection - GE Healthcare's innovative service technology for quick access to technical and clinical experts for personalized support.

InSite with Express Connection enables:

- Proactive monitoring
- Real-time technical and application support**
- Problem diagnosis
- Fast repair of your equipment (more uptime)



1. Please refer to the LOGIQ P9 brochure for a list of purchasable options.

*InSite functions dependent on internet access and speed.

**May vary by region.

LOGIQ and InSite are trademarks of General Electric Company or one of its subsidiaries.



