BRUKSANVISNING


Myrinmätare

OPERATION MANUAL

Goniometer Myrin

KÄYTTÖOHJE

Myrin-mittari

INNEHÅLL

CONTENT

SISÄLTÖ

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Tel. 08-404 12 00

Goniometer Myrin


In order to be able to evaluate the result after treatment of a patient with reduced mobility the passive movement range must be measured before and after the treatment.

The Myrin goniometer makes it possible to measure with precision the movement range in most of the joints of the body. The goniometer consists of a box filled with fluid and the box, which is turnable, is attached to a plate. There is a compass needle in the box that is influenced by the earth magnetism and an inclination needle that is influenced by the gravitation.

The compass needle shows movements on the horizontal plane (around a vertical axis). The inclination needle shows movements on the vertical plane (around a horizontal axis). The box is fastened by means of Velcro-tape to the intended part of the body. Please note that, when moving the goniometer in the room, the compass needle can be affected by iron in the examination table or other objects nearby. During the movement the axis of the goniometer must thus be as close as possible to the movement axis of the joint. This can be made by means of the extension plate.

INSTRUCTIONS OF USE

Fasten the goniometer with the Velcro-tape to the intended part of the body. When measuring around a vertical axis, turn the goniometer in such a way that the compass needle points at zero. When measuring around a horizontal axis, turn the goniometer in such a way that the inclination needle points at zero. Make the movement and read off the scale directly. Check carefully that the needle is moving freely over the whole circuit. The movement is then taking place on the intended level.

EXAMPLES OF APPLICATION

Rotation of the cervical spine

The patient is low-sitting with the Velcro-tape fastened around the head and over the top of it. The goniometer is placed right on the top of the head and the nose pointing forwards. Set the compass needle at zero. Support the back and shoulders of the patient. The patient is actively turning to the right and to the left respectively.

Sideways bending of the cervical spine

Patient with the Velcro-tape fastened around the head. The goniometer is placed at the back of the head or at the front of the head. Set the inclination needle at zero. Patient bending to both sides.

Bending of the cervical spine forwards and backwards

The Velcro-tape is attached around the head. The goniometer is placed on the side of the head. Set the inclination needle at zero. Patient bending the head forwards and backwards respectively.

Rotation of the lumbar back

The patient is sitting astride with the arms in front of the chest and with the hands cupped over the elbows. The goniometer is attached to the angle-bent plate that is placed between the thorax and inner forearm. Set the compass needle at zero. The patient is turning the upper part of the body to the right and to the left respectively.

The above described applications are just examples. Most movements can be measured in different ways, depending on what initial position one chooses. It is in general advantageous to use the inclination needle, i. e. to perform the movement on the vertical plane. You can then disregard of the moving of the goniometer in the room and place it anywhere on the limb. When using the compass needle it must be placed as near as possible to the movement center, as per the explanation in the beginning of this booklet.

SERVICE AND WARRANTY

We offer one year warranty.
For service and spare parts please contact:
Medema Physio AB
E-post: service@medema.se
Tel. +46 8 404 12 00

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Lojer Oy