OXYCYCLE 1

Manual

English.....	2
Nederlands.....	9
Deutsch.....	16
Français.....	23
Espagnol.....	30
Svenska.....	37
Suomi.....	44

IMPORTANT: Read all instructions carefully before using this product. Retain this owner’s manual for future reference. The specifications of this product may vary from this photo and are subject to change without prior notice.

# Table of contents

Important safety precautions .....2

Parts list, hardware, and tools .....3

Exploded view .....4

Assembly instructions.....5

Operating the digital display computer .....6

Operating the tension control knob.....6

Storage precautions .....7

Workout .....7

Warm up and cool down routine.....7

## Important safety precautions


Read all instructions carefully before operating. Retain this owner's manual and keep all original transaction receipts for future reference.

1. Consult your physician before beginning your exercise program.
2. If you experience any irregular physical conditions such as dizziness, severe muscle or joint pain, pain in the chest, stop exercising and consult your physician immediately.
3. Keep children and pets away from the machine while in use.
4. Do not operate the exercise equipment if it is damaged.
5. Set up this equipment as directed in Assembly Instructions and place unit on a flat, stable, and non-slip surface.
6. Keep hands and feet away from any moving parts. Do not insert any object into any openings.
7. Perform 3 to 5 minutes of warm-up and stretching before beginning your exercise.
8. Never exercise to exhaustion.
9. Do not stand on the unit, only exercise in seated position.
10. Keep a minimum safety clearance during operation. The equipment is designed as a home use exercise device. User must have full control of their muscles to operate this unit. **It may not be suitable for people who are quadriplegic or paraplegic.** Please consult your doctor or physician.
11. Close supervision is necessary when this unit is used by, on, or near children, invalids or disabled persons.
12. Use this unit only for its intended use as described in this manual. Do not use attachments not recommended by the manufacturer.
13. Never operate the unit with the air openings blocked; keep the air openings free of lint.
14. Do not use outdoors.
15. Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
16. For Household Use Only.
17. Exercise equipment is to be regularly inspected for maintenance for loose or broken parts.
18. Wear comfortable and suitable clothing when using this equipment. Do not use the equipment barefoot or in socks. Always wear appropriate athletic footwear.
19. Each exercise should be performed in a controlled manner and always begin exercises slowly.
20. Have the unit serviced by a qualified service technician. DO NOT ATTEMPT TO SERVICE THE UNIT YOURSELF.

# Parts list, hardware, and tools

Part #	Description	Q'ty	Part #	Description	Q'ty
001	Main frame	1	026	Hexagon Nut M10x8.5	2
002	Rear Stabilizer Ø50x1.5Tx380	1	027	Tension Bracket	2
003	Front Stabilizer Ø50x1.5Tx220	1	028	Hexagon Nut M6x6	2
004	Magnet Bracket	1	029	Hexagon Nut M8x7.7	1
005	Tension Control Knob Bracket	1	030	Washer Ø8.5x16x1.5T	1
006	Belt Pulley Ø120x152	1	031	Hexagon Socket Button Head Bolt M8x42	1
007	Magnet 15x10x5T	7	032	Cross Recessed Truss Head Tapping Screw ST4.2x15	4
008	Spacer Ø14xØ10.5x4	2	033	Cross Recessed Pan Head Tapping Screw ST4.2x15	7
009	Flywheel	1	034	Single Sided Adhesive Tape	4
010	Stabilizer End Cap	4	035	PVC Pad Ø18.5xØ17x9	2
011	Square Sensor (L=300mm)	1	036	Cross Recessed Pan Head Tapping Screw ST2.9x10	2
012	Big Washer Ø75x2.0T	2	037	Spring Ø9.4xØ10x50	1
013	Bearing Bracket Ø74xØ41.5x12.5	2	038	Belt 220PJ3	1
014	Right Crank 5" (9/16")	1	039	Crank Cap	2
015	Left Crank 5" (9/16")	1	040	Right Cover	1
016	Digital Display Computer	1	041	Left Cover	1
017	Hexagon Flange Bolt M8x20	2	042	Right Hand grip 9/16"	1
018	Hexagon Socket Button Head Bolt M8x12	6	043	Left Hand grip 9/16"	1
019	Spring Washer Ø8.4xØ13.5x2.5T	6	044	Right Foot Pedal (JD-5B)	1
020	C Ring Ø17xØ15.7x1.0T	4	045	Left Foot Pedal (JD-5B)	1
021	Bearing 60O3Z	2	046	Round Magnet Ø12x5.0t	1
022	Tension Control Knob	1	047	Non-Slip Mat	1
023	Flywheel Axle Ø10x75	1	048	Cross Recessed Truss Head Tapping Screw ST4.2x20	2
024	Cross Recessed Truss Head Bolt M5x12	10	049	Hexagon Nut M5x4	1
025	Eyebolt M6x43	2	050	Cross Recessed Pan Head Bolt M5x15	1

(18) Hexagon Socket  
Button Head Bolt  
6 PCS

(19) Spring Washer  
6 PCS


Wrench  
1 PC




Allen Wrench  
1 PC

# Exploded view

## Assembly instructions

1. Hold the unit in an upright position. The tension control knob (22) side facing up. Then insert the front stabilizer (3) into the main frame (1).
2. Attach the front stabilizer (3) to the main frame (1) with three hexagon socket button head bolts (18) and three spring washers (19). Tighten bolts with the allen wrench provided. (p.3)
3. Insert the rear stabilizer (2) into the main frame (1).
4. Attach the Rear Stabilizer (2) to the Main Frame (1) with three Hexagon Socket Button Head Bolts (18) and three Spring Washers (19). Tighten bolts with the Allen Wrench provided. (p.3)
5. Connect the left foot pedal (45) to the left crank (15) assembly where you would position your left foot when exercising. Thread it into the crank assembly in a counter-clockwise direction with the wrench provided. (See figures 5 and 6).
6. Connect the right foot pedal (44) to the right crank (14) assembly by threading it in a clockwise direction.

**Note:** The left/right pedals and left/right cranks are marked with "L" and "R" (left and right).

## Operating the digital display computer

### Computer Button Functions

Push the button to turn on the display.

Push the button to select a function.

Push and hold the button for 3 seconds to reset all functions to zero except for the TOTAL REPS function.


### Computer Display Functions & Measuring Range

SCAN	Automatically displays each function every 6 seconds.
SCAN → REP → TMR → T.REPS → CAL → DIS	
REPS	Displays the number of time of rotation.
TMR	Displays the exercise time.
T.REPS	Displays the number of total rotation.
CAL	Displays the calories burned during exercise.
DIS	Displays the accumulative distance traveled during exercise.
STOP	The unit has stopped.
AUTO START	Start motion or push button.
AUTO DISPLAY SHUT OFF	After approximately 4 minutes without operation.

### Battery replacement

Remove the display by using a slotted screwdriver.

Remove the battery with a slotted screwdriver. Replace with a new 3V lithium button cell battery.

## Operating the tension control knob

### Tension control knob

The resistance can be adjusted by turning the resistance control knob. To increase the resistance, turn the resistance control knob in a clockwise direction. To decrease the resistance, turn the resistance control knob in a counterclockwise direction.

### How to use the OxyCycle 1 Pedal Exerciser

When using the OxyCycle 1, the proper way to exercise is to set the desired resistance by turning the tension control knob, and then pedaling to get the best results. Pedaling at excessive speeds or applying excessive force will damage the product or may cause injury to user. Please start exercising at a slow pedaling speed.

## Storage precautions

Do not store the unit in any place where it will be subjected to high or low temperature.  
Do not expose to direct sunlight or outdoor for extended periods of time.  
Do not store the unit in a humid or dusty area.

## Workout

### **Leg exercise**

Place the Pedal Exerciser on a non-slip mat and sit comfortably in a chair, the device directly in front of you. The knee joint should not be flexed less than 90 degrees. Wear rubber soled shoes and adjust foot pedal straps to the proper length. Barefooted use of the OxyCycle is not recommended.

### **WARNING**


Do not stand on the unit to exercise. Only exercise in a seated position. Make sure the seat does not move or pivot. Do not sit in a chair with casters.

NOTE: Use the foot pedals for leg exercises.

### **Arm Exercise**

Place the unit directly in front of you on a table top (on the included anti-slip mat). Use the hand grips for arm exercises.

## Warm up and cool down routine

The WARM-UP is an important part of any workout. The purpose of warming up is to prepare your body for exercise and to minimize injuries. Warm up for two to five minutes before aerobic exercising. It should begin every session to prepare your body for more strenuous exercise by heating up and stretching your muscles, increasing your circulation and pulse rate, and delivering more oxygen to your muscles.

COOL DOWN at the end of your workout, repeat these exercises to reduce soreness in tired muscles. The purpose of cooling down is to return the body to its resting state at the end of each exercise session. A proper cool-down slowly lowers your heart rate and allows blood to return to the heart.

### 1. Head rolls

Rotate your head to the right for one count, you should feel a stretching sensation up the left side of your neck. Then rotate your head back for one count, stretching your chin to the ceiling and letting your mouth open. Rotate your head to the left for one count, then drop your head to your chest for one count.

### 2. Shoulder lifts

Lift your right shoulder toward your ear for one count. Then lift your left shoulder up for one count as you lower your right shoulder.

### 3. Side stretches

Open your arms to the side and lift them until they are over your head. Reach your right arm as far toward the ceiling as you can for one count. Repeat this action with your left arm.

### 4. Quadriceps stretch

With one hand against a wall for balance, reach behind you and pull your right foot up. Bring your heel as close to your buttocks as possible. Hold for 15 counts and repeat with left foot.

### 5. Inner thigh stretch

Sit with the soles of your feet together and your knees pointing outward. Pull your feet as close to your groin as possible. Gently push your knees toward the floor. Hold for 15 counts.

### 6. Toe touches

Slowly bend forward from your waist, letting your back and shoulders relax as you stretch toward your toes. Reach as far as you can and hold for 15 counts.

### 7. Hamstring stretches

Extend your right leg. Rest the sole of your left foot against your right inner thigh. Stretch toward your toe as far as possible. Hold for 15 counts. Relax and then repeat with left leg.

### 8. Calf/achilles stretch

Lean against a wall with your right leg in front of the left and your arms forward. Keep your left leg straight and the right foot on the floor; then bend the right leg and lean forward by moving your hips toward the wall. Hold, then repeat on the other side for 15 counts.

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Qarad UK Ltd.