





MAIN FEATURES

- Manual power control
- No Power adapter is needed

THE BASICS

- Good cycling feeling
- 5 kHz HR system (only in display)
- High quality components for heavy duty users
- Suitable for training and rehabititation

Max user weight 250 kg

Medical Device class IIa



Specifications MONARK 827E

Certifications	FU directive
Sor timoutions	93/42/EU class IIa
Frame	
Туре	Monark Classic
Handlebar	Monark handlebar
Stem	Monark Standard
Seatpost	Monark Standard
Saddle	Saddle, blue
ВВ	68/122 mm
Crank	Steel, 170 mm
Pedals	9/16", pedals with foot straps
RESISTANCE	
System	Monark Basic
Flywheel	20 kg
DISPLAY	
RPM	Yes
HR	5 kHz
Time	Yes
Speed	Yes
Distance	Yes
Watt	n/a
Connectivity	
RS232	n/a
Usb	n/a
WEIGHT & MEASUREMENTS	
Max user weight	250 kg
Lenght (mm)	1120 mm
Width (mm)	519 mm at handlebar 640 mm at support legs
Height (mm)	945-1295 mm at handlebar 780-1105 mm at saddle
Weight (kg)	47 kg
POWER REQUIREMENTS	Batteries in the display



HEART RATE

Equipped to measure heartrate via chestbelt or sensors in the handles.

MEDICAL CE

Monark Exercise AB is since 1992 ISO 9001 certified. All models are:

- Approved under CE Certificate No. 153 601 (MDD IIa)
- Anti-corrosion treated and powder coated
- Made in Sweden

Upplaga 1. 2016-02-28

