



MAIN FEATURES

- Compact design
- Genuine cycling feeling
- No Power adapter needed

THE BASICS

- User friendly display
- 5 kHz HR system
- Watt in display

Max user weight 250 kg

Medical Device class IIa



Specifications MONARK 915E

Certifications	EU directive 93/42/EU class IIa
Frame	
Туре	Monark Classic Mini
Handlebar	Monark handlebar
Stem	n/a
Seatpost	Monark Standard
Saddle	Saddle Blue
ВВ	68/122 mm
Crank	Steel, 170 mm
Pedals	9/16", pedals with foot straps
RESISTANCE	
System	Monark Basic
Flywheel	17 kg
DISPLAY	
RPM	n/a
HR	5 kHz
Time	Yes
Speed	Yes
Distance	Yes
Watt	Yes
Connectivity	
RS232	n/a
Usb	n/a
WEIGHT & MEASUREMENTS	
Max user weight	125 kg
Lenght (mm)	870 mm
Width (mm)	530 mm at handlebar 460 at support legs
Height (mm)	800 mm at seat
Weight (kg)	38 kg
POWER REQUIREMENTS	Batteries in the display



HEART RATE

Equipped to measure heartrate via chestbelt or sensors in the handles.

MEDICAL CE

Monark Exercise AB is since 1992 ISO 9001 certified. All models are:

- Approved under CE Certificate No. 153 601 (MDD IIa)
- Anti-corrosion treated and powder coated
- Made in Sweden

Upplaga 1. 2016-02-28

