



MONARK

828E



MAIN FEATURES

- The most sold Ergometer cycle ever now with novo mini display
- Mechanic calibration
- Manual power control
- Power bank

THE BASICS

- Good cycling feeling
- Easy to calibrate
- ANT+ HR system
- ANT+ with FE-C protocol (Fitness Equipment Controls)
- BLE with FTMS protocol (Fitness Machine Service)

Max user weight 250 kg

Medical Device class IIa

Interactive cycling

Specifications

MONARK 828E

Certifications	EU directive 93/42/EU class IIa
Frame	
Type	Monark Classic
Handlebar	Monark Handlebar
Stem	Monark Standard
Seatpost	Monark Standard
Saddle	Saddle, blue
BB	68/122 mm
Crank	Steel, 170 mm
Pedals	9/16", pedals with foot straps
RESISTANCE	
System	Pendulum
Flywheel	20 kg
DISPLAY	
RPM	Yes
HR	ANT+
Time	Yes
Speed	Yes
Distance	Yes
Watt	Yes
Connectivity	
ANT+	FE-C (Fitness Equipment Controls)
BLE	FTMS (Fitness Machine Service)
WEIGHT & MEASUREMENTS	
Max user weight	250 kg
Length (mm)	1120 mm
Width (mm)	518 mm at handlebar 640 mm at support legs
Height (mm)	945-1295 mm at handlebar 780-1105 mm at saddle
Weight (kg)	52 kg
POWER REQUIREMENTS	Power bank



HEART RATE

Equipped to measure heartrate via chestbelt.



TEST

Offers the possibility of performing different tests thanks to the accurate workload which can be calibrated.

MEDICAL CE

Monark Exercise AB is since 1992 ISO 9001 certified. All models are:

- Approved under CE Certificate No. 15 36 01 (MDD IIa)
- Anti-corrosion treated and powder coated
- Made in Sweden